

NEW AlteCap™ Solo

Single screw cap handling for any manufacturer's tube types

Quickly cap, or recap single tubes. Compatible with AlteTubes and other manufacturer's screw cap tubes in SBS rack format.

Part No. Product name

404001 AlteCap Solo



Flexible

Easy to use, intuitive touch button controls. Light and portable

Precise torque control


Internal cap driver control ensures each screw cap is perfectly capped

Functionality

Recap tubes or discard used caps. Apply fresh caps from a cap carrier

Multi use stand

Convenient for storage, USB charging and holding caps for re-capping



AlteCap™

Screw cap handling

Flexible systems for screw capping tasks

Automating screw cap decapping and capping is now easier and more affordable with AltemisLab's comprehensive range of decapping and capping systems. Accommodating all throughput requirements, the AlteCap™ decappers can help with the time challenges screw cap tube users often experience. All systems in the AlteCap™ range are compatible with a variety of different tube types and rack formats from any manufacturer.

Transforming lab technique.


AlteCap™ Solo



AlteCap™ Switch



AlteCap™ Swift



Tube throughput	1	6, 8 & 12	24, 48 or 96
Operation modes	Cap – Decap – Recap	Cap – Decap – Recap	Cap – Decap – Recap
Type	Handheld	Handheld	Benchtop
Compatible with SBS racks	Yes	24, 48 or 96	24, 48 or 96
Tube manufacturer	Any	Any	Any
Fresh capping	Cap carrier	Cap carrier	Cap carrier
Stand included	Yes	Yes	N/A
Integration friendly	N/A	N/A	Yes
Dimensions - cm	2.5 x 19	19 x 12.5 x 82	36.5 x 24.7 x 33.2
Power	USB charging	USB charging	Mains powered